

POLYESTER weiß

TIPOLOGY:	Braided non-absorbable synthetic suture.
COLORS:	White.
MOLECULE:	Polyester.
CALIPER:	USP from 12/0 to 8; EP from 0,1 to 9.
COATING TYPE:	Uncoated.
SPECIFICS:	<ul style="list-style-type: none">• Constant tensile strength through time.• Easy handling, flexibility and easy execution of the node.• Excellent knot tying security.• Maximum reliability in the application of prosthesis and in the approximation of tissue subjected to tension.• Reduced capillarity.
INFORMATIONS:	For its particular characteristics, Polyester weiß sutures are especially suitable for cardiovascular surgery and suturing of the skin. It can also be employed in the approximation and/or ligation of the soft tissues, in ophthalmic and orthopedic surgery, and for neural tissues.
SUTURE TENSILE STRENGTH:	Polyester weiß doesn't show significant loss of strength through time.
ABSORPTION:	Non-absorbable.
STERILIZATION:	Ethylene oxide.
EXPIRE:	5 years.
CERTIFICATION:	<ul style="list-style-type: none">- CE 0373 Notified Body "Istituto Superiore di Sanità" per class III.- EPG-0164-18 Addendum n.01-19.- EPG-0164-18 Addendum n.02-20.- QCT-0087-18 Addendum n.03-19.- QCT-0087-18 Addendum n.04-20.
MANUFACTURER:	BIO SUD MEDICAL SYSTEMS S.R.L. Via dei Fabbri 23-25, 70026 - Modugno (BA) - Italy